

Innova Pharm

FUNCTIONAL EXPERTISE SINCE 1999



MADE IN FRANCE

Our devices are manufactured in France using medical-grade materials that meet international standards (FDA, USP Class VI), for safe, controlled use.

Formulated without bisphenols or phthalates.

Over

28 years

of expertise

Over

80 devices

for orofacial myofunctional therapy

Up to

5 densities

for optimal patient comfort

Advice and indications



Adapting to the device

Sensitivity in the front teeth may occur on waking. It is temporary, subsides quickly and reflects the device's action on the incisors and canines. During the first few weeks (up to around 4 weeks), it is common for the appliance to fall out during the night. This is normal: functional retraining takes time and perseverance. The devices are not recommended for patients with sleep apnoea.

Cleaning the devices



- ▶ Clean the device before placing it in the mouth and after each removal, using a dedicated soft toothbrush and a mild cleaning agent (toothpaste, neutral soap, etc.).
- ▶ Gently brush all surfaces, then rinse thoroughly with clean water.
- ▶ In addition, once a week, soak the device for 30 minutes in a container of water mixed with mouthwash or an equivalent solution.
- ▶ After soaking, brush again, then rinse thoroughly with clean water before use.
- ▶ Dry the device thoroughly before storing it in its case, making sure the case remains clean.

Sterilisation



The Écran and Déglutix can be sterilised in an autoclave at 134 °C using a PRION cycle.

Never boil the device.

Hygiene and precautions



Poor oral hygiene or inadequate cleaning of the device and its case may cause inflammatory reactions. It is essential to maintain strict hygiene throughout treatment.

A range designed to support orofacial myofunctional therapy

FIRST-LINE AND EXERCISE DEVICES	8
Breathing and swallowing retraining	
SCREEN	10
DÉGLUTIX	12
INTERCEPTIVE DEVICES	14
Correcting orofacial myofunctional disorders	
CLASS I	18
CLASS II	20
CFP APPLIANCE	22
SQUARE APPLIANCE	23
CLASS III	24
AUXILIARY TREATMENT DEVICES	26
Supporting orthodontic treatment	
ALIGNER APPLIANCE	28
BRACES APPLIANCE	29
RETENTION DEVICES	30
Maintaining results	
PROCCLUSEUR	32

CONTENTS

Key



Recommended minimum age



Clinical observation



Compatible with orthodontic treatment



Daytime wear



Night-time wear

Prices⁽¹⁾

- ✓ Volume discounts & mix-and-match all year round ⁽²⁾
- ✓ 1 FREE device for every 10 ordered ⁽³⁾
- ✓ Free shipping on orders of 10 or more devices⁽⁴⁾

FIRST-LINE AND EXERCISE DEVICES

SCREEN ⁽⁵⁾					
Size	Density		x5	x20	Page
3-6 yrs			13 €	12 €	10
6+ years			13 €	12 €	10
6+ years	RIGID		13 €	12 €	10
3-6 yrs		WITH BITE STOP	14 €	13 €	11
6+ years		WITH BITE STOP	14 €	13 €	11

DÉGLUTIX					
Size	Colour	Unit	x5	x20	Page
3-6 yrs	Blue	35 €	25 €	22 €	12
3-6 yrs	Pink	35 €	25 €	22 €	12
3-6 yrs	Orange	35 €	25 €	22 €	12
6+ years	Blue	35 €	25 €	22 €	12
6+ years	Pink	35 €	25 €	22 €	12
6+ years	Orange	35 €	25 €	22 €	12
6+ years	Blue	WITH RAMPS	35 €	25 €	13
6+ years	Pink	WITH RAMPS	35 €	25 €	13
6+ years	Orange	WITH RAMPS	35 €	25 €	13
XL	Blue		35 €	25 €	12
XL	Pink		35 €	25 €	12
XL	Orange		35 €	25 €	12

INTERCEPTIVE DEVICES

CLASS I							
Size	Density	Unit	x5	x10	x20	x50	x100 Page
S	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 18-19
S	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 18-19
S	RIGID	39 €	35 €	30 €	28 €	25 €	23 € 18-19
M	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 18-19
M	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 18-19
M	RIGID	39 €	35 €	30 €	28 €	25 €	23 € 18-19
M	SOFT MEDIUM	46 €	42 €	39 €	35 €	33 €	31 € 18-19
M	MEDIUM RIGID	46 €	42 €	39 €	35 €	33 €	31 € 18-19
L	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 18-19
L	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 18-19
L	RIGID	39 €	35 €	30 €	28 €	25 €	23 € 18-19
XL	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 18-19
XL	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 18-19

CLASS II							
Size	Density	Unit	x5	x10	x20	x50	x100 Page
S	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 20-21
S	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 20-21
S	RIGID	39 €	35 €	30 €	28 €	25 €	23 € 20-21
S	SOFT MEDIUM	46 €	42 €	39 €	35 €	33 €	31 € 20-21
S	MEDIUM RIGID	46 €	42 €	39 €	35 €	33 €	31 € 20-21
M	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 20-21
M	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 20-21
M	RIGID	39 €	35 €	30 €	28 €	25 €	23 € 20-21
M	SOFT MEDIUM	46 €	42 €	39 €	35 €	33 €	31 € 20-21
M	MEDIUM RIGID	46 €	42 €	39 €	35 €	33 €	31 € 20-21
L	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 20-21
L	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 20-21
L	RIGID	39 €	35 €	30 €	28 €	25 €	23 € 20-21
L	SOFT MEDIUM	46 €	42 €	39 €	35 €	33 €	31 € 20-21
L	MEDIUM RIGID	46 €	42 €	39 €	35 €	33 €	31 € 20-21
XL	SOFT	39 €	35 €	30 €	28 €	25 €	23 € 20-21
XL	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 € 20-21

(1) Prices are in euros, inclusive of tax, subject to changes in applicable taxes (including VAT) and/or selling prices. Prices are valid from 1 October 2026 and may change during the year due to changes in applicable taxes or economic conditions. Before ordering, check the latest prices at www.innovapharm.fr.
(2) Mix-and-match orders are available for devices with the same price.
(3) One free device for every 10 ordered.
(4) Free shipping within mainland France for orders of 10 or more devices. Offer applies to devices only and cannot be combined with other offers.
(5) Ecran devices are sold in quantities of at least 5.

CFP APPLIANCE

Size	Density	Unit	x5	x10	x20	x50	x100	Page
ONE SIZE	SOFT	39 €	35 €	30 €	28 €	25 €	23 €	22
ONE SIZE	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	22
ONE SIZE	RIGID	39 €	35 €	30 €	28 €	25 €	23 €	22

SQUARE APPLIANCE

Size	Density	Unit	x5	x10	x20	x50	x100	Page
N°1	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°2	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°3	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°4	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°5	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°6	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°7	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°8	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°9	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23
N°10	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	23

CLASS III

Size	Density	Unit	x5	x10	x20	x50	x100	Page
S	SOFT	39 €	35 €	30 €	28 €	25 €	23 €	24-25
S	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	24-25
M	SOFT	39 €	35 €	30 €	28 €	25 €	23 €	24-25
M	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	24-25
M	SOFT MEDIUM	46 €	42 €	39 €	35 €	33 €	31 €	24-25
M	MEDIUM RIGID	46 €	42 €	39 €	35 €	33 €	31 €	24-25
L	SOFT	39 €	35 €	30 €	28 €	25 €	23 €	24-25
L	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	24-25

AUXILIARY DEVICES FOR TREATMENT

ALIGNER APPLIANCE

Size	Density	Unit	x5	x10	x20	x50	x100	Page
M	MEDIUM	39 €	35 €	30 €	28 €	25 €	23 €	28

BRACES APPLIANCE

Size	Density		Unit	x5	x10	x20	x50	x100	Page
M	MEDIUM		39 €	35 €	30 €	28 €	25 €	23 €	29
M	MEDIUM	WITH BLOCKS	39 €	35 €	30 €	28 €	25 €	23 €	29
M	MEDIUM	SINGLE ARCH	39 €	35 €	30 €	28 €	25 €	23 €	29

RETENTION DEVICES

PROCCLUSEUR

Size		Unit	x5	x10	x20	x50	x100	Page
N°1		39 €	35 €	30 €	28 €	25 €	23 €	32
N°2		39 €	35 €	30 €	28 €	25 €	23 €	32
N°3		39 €	35 €	30 €	28 €	25 €	23 €	32
N°4		39 €	35 €	30 €	28 €	25 €	23 €	32
N°5		39 €	35 €	30 €	28 €	25 €	23 €	32
N°6		39 €	35 €	30 €	28 €	25 €	23 €	32
N°7		39 €	35 €	30 €	28 €	25 €	23 €	32
N°8		39 €	35 €	30 €	28 €	25 €	23 €	32
N°1	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33
N°2	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33
N°3	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33
N°4	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33
N°5	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33
N°6	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33
N°7	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33
N°8	ADJUSTABLE	46 €	42 €	39 €	35 €	33 €	31 €	33

Sizes and densities

A prescription tailored to each patient

Innova Pharm devices are available in several sizes and densities to match the choice of device to the patient's functional needs, comfort and tolerance.

Choosing the right size: The size should ensure stable, comfortable positioning, effective functional action and good patient compliance.

FIRST-LINE & EXERCISE

3-6
years

6+
years

XL

Adolescents and
adults



INTERCEPTIVE DEVICES

S

3-6 years

M

6+ years

L

Includes second
molars

XL

Adolescents and
adults



SQUARE APPLIANCE

Size is determined by adding the mesiodistal widths of the 4 upper incisors and using a size conversion chart (1).



Scan this QR
code for the
step-by-step
method.



PROCCLUSEUR⁽²⁾

Learn how to use the measuring ruler (3) and select the most suitable size for your patient.



Scan the QR
code.



(1) Size conversion chart on page 23. Also available on request when purchasing a Square Appliance. (2) More information on Proccluseur and adjustable Proccluseur devices on pages 32-33. (3) Measuring ruler available on request when purchasing a Proccluseur or adjustable Proccluseur.

For more information on choosing sizes,
scan this QR code to watch our
explanatory video.



Understanding densities

Densities correspond to different levels of material hardness. They allow comfort, muscle control, therapeutic action and functional stability to be adjusted.

SOFT DENSITY

Comfort & acceptance

Young patients, hypersensitivity;
recommended for size S

MEDIUM DENSITY

Therapeutic versatility

Most treatments; recommended
for size M and above

RIGID DENSITY

Enhanced control

High muscle tone; need for
greater control

Innova Pharm dual-density devices

Patented technology combining two density levels in one device: softness in contact with soft tissues and rigidity to control the dental arches, without increasing thickness.



Prescribing guidance

COMFORT FIRST

SOFT DENSITY

VERSATILE USE

MEDIUM DENSITY

ENHANCED CONTROL

RIGID DENSITY

COMFORT + EFFECTIVENESS

DUAL DENSITY SOFT MEDIUM

CONTROL + STABILITY

DUAL DENSITY MEDIUM RIGID



The device, its size and its density must be prescribed by a qualified healthcare professional. Further information on sizes and densities is provided on the dedicated product pages.

For more information on choosing our densities, scan this QR code to watch our explanatory video.



First-line and exercise devices

Breathing and swallowing retraining

These orofacial myofunctional devices provide a first therapeutic step in managing orofacial dysfunctions in children.




SCREEN

Breathing retraining

Designed for breathing retraining, the Écran helps children develop exclusively nasal breathing to support harmonious development of the upper and lower jaws.



Developing nasal breathing



The Écran sits between the lips and teeth and applies pressure to the upper incisors when there is an anteroposterior discrepancy between the teeth. It also reduces cheek contraction, which causes narrow palates.

It provides an alternative to finger sucking and pacifier use. It can help eliminate this harmful habit and restore a functional environment that supports harmonious growth of the jaws and teeth.

Available in:

3-6
years

6+
years

RIGID DENSITY

Recommended for patients with high lip and cheek muscle tone.

Available in:

6+
years



i THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

SCREEN

WITH BITE STOP

| Breathing retraining

Recommended for dental discrepancies where the lower jaw is set back.



Bite stop

Encourages mandibular advancement and directs the tongue towards the palatal papillae for mature swallowing. The bite stop sits behind the lower incisors.

Available in:

3-6
years6+
years

Daytime wear gradually prepares the patient for night-time wear, improving acceptance and reducing the likelihood of the device falling out during the night.



FUNCTIONAL RETRAINING COMBINED WITH ORTHODONTIC TREATMENT CAN SUPPORT CLINICAL PROGRESS.

CLINICAL OBSERVATION

Developing nasal breathing

Clinical case showing the progressive evolution of an anterior open bite during orofacial myofunctional retraining with the Écran.



APRIL 1998



JANUARY 2000

 Prescription and follow-up must be provided by a qualified healthcare professional.



Not suitable for children under 36 months, patients with Class III occlusion or prognathic tendencies, or exclusive mouth breathers.


DÉGLUTIX

| Swallowing retraining

Designed to retrain immature swallowing through reflex action, Déglutix promotes correct, natural swallowing without conscious effort.



Developing mature swallowing



Déglutix sits between the lips and teeth. To swallow, the patient brings the teeth together and places the tongue against the palate to achieve correct mature swallowing, assisted by the air vents. This new swallowing pattern occurs without conscious effort, allowing learning to take place spontaneously.

Available in:

3-6
years

6+
years

XL

Colours available by size:

Blue

Orange

Pink



i THE DEVICE MUST BE WORN FOR APPROXIMATELY 30 MINUTES A DAY.

DÉGLUTIX

WITH RAMPS

| Swallowing retraining

Recommended
for occlusal
instability.

Available in:



6+
years

Available colours:

Blue

Orange

Pink



FUNCTIONAL RETRAINING COMBINED WITH ORTHODONTIC TREATMENT CAN SUPPORT CLINICAL PROGRESS.

CLINICAL OBSERVATION

Results in less than 3 months

In a patient receiving fixed-brace treatment, wearing Déglutix alongside their care led to functional improvement associated with a change in the open bite.



DECEMBER 2016



FEBRUARY 2017

 Prescription and follow-up must be provided by a qualified healthcare professional.



Not suitable for children under 36 months.

Interceptive devices

Correcting orofacial **myofunctional disorders**

These devices are used as a second-line treatment.
They correct functional and occlusal imbalances
while supporting harmonious development of the dental
arches and underlying bone structures.




Our interceptive devices act simultaneously on breathing,
tongue posture, occlusion and muscle balance,
for comprehensive functional care.

Vestibular shield

Corrects mouth breathing.

Lip bumper extensions

Neutralise lip and chin contractions, reduce muscle pressure and promote mandibular growth.



Patented bilateral blocks

Innova Pharm blocks act as a tongue elevator. They position the tongue against the palate and promote mature swallowing.

U-shaped lower arch

Fits within the perimeter of the upper arch to encourage interdigitation of the posterior teeth.

Occlusal plane

Guides the relationship between the arches and contributes to occlusal correction.

Our densities

To tailor each device to the patient's functional needs and comfort, we offer several density levels without increasing thickness:

SOFT DENSITY

MEDIUM DENSITY

RIGID DENSITY

Each density addresses a specific indication, allowing therapeutic action to be adjusted progressively while supporting patient compliance.

Unlike conventional approaches that add material to increase rigidity, our devices maintain a consistently slim profile for greater comfort and an optimal fit in the mouth.

DUAL DENSITY: EXCLUSIVELY BY INNOVA PHARM

Available in:

DUAL DENSITY SOFT MEDIUM

DUAL DENSITY MEDIUM RIGID

Comfort and effectiveness combined

The first appliances to combine two densities in a single device, using an innovative, patented manufacturing process.

Combines:

- Softness against soft tissues (lips, cheeks)
- Rigidity at the dental arches

For greater patient comfort and optimised therapeutic effectiveness.

Designed for clinical practice

FOR THE CLINICIAN

Precise adaptation to each indication

Controlled treatment progression

Optimised results

FOR THE PATIENT

Greater everyday comfort

Better compliance

Easier acceptance

CLASS I

Maintenance and balance in Class I

For Class I occlusion, this appliance helps correct tongue and lip dysfunctions.



- Vestibular shield ✓
- Patented bilateral blocks ✓
- Lip bumper extensions ✓
- Occlusal plane ✓
- U-shaped arch ✓

Key feature

The lower jaw remains in a Class I position.

Orofacial myofunctional and occlusal retraining.

Main actions

- Corrects tongue and lip dysfunctions
- Improves nasal breathing
- Retrains swallowing
- Stabilises occlusal balance
- Relieves temporomandibular joints (TMJs)

Instructions for use

The bilateral blocks point downwards beneath the tongue. The tongue naturally adopts a raised position.

THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

Maintenance and balance in Class I



Available in:

4 sizes - 5 densities

	S	M	L	XL
SOFT DENSITY	✓	✓	✓	✓
MEDIUM DENSITY	✓	✓	✓	✓
RIGID DENSITY	✓	✓	✓	
DUAL DENSITY SOFT MEDIUM		✓		
DUAL DENSITY MEDIUM RIGID		✓		

CLINICAL OBSERVATION

In 8-year-old twins of comparable age and environment, differences in breathing function can have visible effects on orofacial balance and dental arch development.



Twin A



Twin B

Training in functional therapy*

Explore our courses on managing orofacial dysfunctions.



Scan this QR code to register for our courses.

* More information on our training courses on page 35.

CLASS II

| Correcting Class II discrepancies

For Class II occlusion, this appliance helps correct tongue and lip dysfunctions while encouraging mandibular advancement.



- Vestibular shield ✓
- Patented bilateral blocks ✓
- Lip bumper extensions ✓
- Occlusal plane ✓
- U-shaped arch ✓

Key feature

The lower jaw is positioned edge to edge to encourage interdigitation of the posterior teeth.

Corrects anteroposterior discrepancies.

Main actions

- Corrects tongue and lip dysfunctions
- Improves nasal breathing
- Retrains swallowing
- Stabilises occlusal balance
- Relieves temporomandibular joints (TMJs)

Instructions for use

The bilateral blocks point downwards beneath the tongue. The tongue naturally adopts a raised position.

THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

Correcting Class II discrepancies



Available in:

4 sizes - 5 densities

	S	M	L	XL
SOFT DENSITY	✓	✓	✓	✓
MEDIUM DENSITY	✓	✓	✓	✓
RIGID DENSITY	✓	✓	✓	
DUAL DENSITY SOFT MEDIUM	✓	✓	✓	
DUAL DENSITY MEDIUM RIGID	✓	✓	✓	

CLINICAL OBSERVATION

In a patient with Class II occlusion, use of the appliance improved function and progressively corrected the discrepancy between the upper and lower jaws.



FEBRUARY 2022



OCTOBER 2022

Training in functional therapy*

Explore our courses on managing orofacial dysfunctions.



Scan this QR code to register for our courses.

* More information on our training courses on page 35.

CFP APPLIANCE

Correcting Class II discrepancies

A simple, immediate solution for managing dysfunctions in Class II occlusion.



Vestibular shield



Patented bilateral blocks



Occlusal plane



U-shaped arch



Available in one size
and 3 densities*:

One
size

SOFT DENSITY

MEDIUM DENSITY

RIGID DENSITY

Key feature

The lower jaw is positioned edge to edge.

The first Class II appliance designed by Innova Pharm, the CFP is an intermediate size between S and M*.

Main actions

- Corrects tongue and lip dysfunctions
- Improves nasal breathing
- Retrains swallowing
- Stabilises occlusal balance
- Relieves temporomandibular joints (TMJs)

Instructions for use

The bilateral blocks point downwards beneath the tongue. The tongue naturally adopts a raised position.

THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

Prescription and follow-up must be provided by a qualified healthcare professional.



Not suitable for children under 36 months or patients with sleep apnoea.

* The CFP appliance uses a specific formulation that differs from other Innova Pharm devices, while retaining comparable indications and instructions for use.

SQUARE ARCH APPLIANCE

PLANAS tracks
Correcting Class II
discrepancies

Designed according to the KORKHAUS chart for use after PLANAS track treatment.



Vestibular shield



Patented bilateral blocks



Occlusal plane



Square arch shape



Lip bumper extensions



How to measure?

Determine the size using the matching table below, based on the KORKHAUS chart. If the measurement falls between two sizes, choose the larger size.



Scan this QR code to watch the step-by-step explanatory video.

Appliance size	1	2	3	4	5	6	7	8	9	10
Sum of mesiodistal widths of the 4 upper incisors (SI)	25	26.2	27.5	28.7	30	31.2	32.5	33.7	35	36.2

Key feature

Its square upper arch shape, designed according to the KORKHAUS chart, promotes posterior occlusion. The advancement plane frees the lower jaw and encourages mature swallowing. The lower jaw is positioned edge to edge.

Main actions

- Corrects tongue and lip dysfunctions
- Improves nasal breathing
- Retrains swallowing
- Stabilises posterior occlusal balance
- Relieves temporomandibular joints (TMJs)
- Placement and follow-up of PLANAS tracks

Instructions for use

Place the appliance in the mouth, aligned with the upper arch. Determine its size using a matching table* based on the KORKHAUS chart. Wearing time is set by the practitioner according to the treatment protocol.

- 1 THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.**



Prescription and follow-up must be provided by a qualified healthcare professional.



Not suitable for children under 36 months or patients with sleep apnoea.

* A matching table is available on request when purchasing the Square Arch Appliance.

CLASS III

Correcting Class III discrepancies
Anterior crossbites

For non-skeletal Class III occlusion or an anterior crossbite, this appliance helps correct orofacial dysfunctions while improving the relationship between the incisors.



- Vestibular shield ✓
- Patented bilateral blocks ✓
- Lip bumper extensions ✓
- Occlusal plane ✓
- U-shaped arch ✓

Key feature

Incisal overjet is increased to correct anterior crossbites.

The appliance corrects Class III discrepancies of functional origin.

Main actions

- Corrects tongue and lip dysfunctions
- Improves nasal breathing
- Retrains swallowing
- Stabilises occlusal balance
- Relieves temporomandibular joints (TMJs)
- Corrects negative incisal overjet

Instructions for use

The bilateral blocks point downwards beneath the tongue. The tongue naturally adopts a raised position.

THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

Corrects anterior crossbite



Available in:

3 sizes - 4 densities

	S	M	L
SOFT DENSITY	✓	✓	✓
MEDIUM DENSITY	✓	✓	✓
DUAL DENSITY SOFT MEDIUM		✓	
DUAL DENSITY MEDIUM RIGID		✓	

CLINICAL OBSERVATION

In a patient with a crossbite, wearing the Class III appliance progressively corrects incisal overjet and improves occlusal relationships.



DECEMBER 2014



JULY 2015

Training in functional therapy*

Explore our courses on managing orofacial dysfunctions.



Scan this QR code to register for our courses.

* More information on our training courses on page 35.

Adjunctive treatment appliances

Oro-myofunctional support
during treatment with
fixed braces and aligners.

These appliances complement treatment to optimise its effectiveness and retrain orofacial functions.



ALIGNER APPLIANCE

Adjunctive appliance for
aligner treatment

Designed to integrate oromyofunctional retraining into
orthodontic aligner treatment.



Available in:

M

MEDIUM DENSITY

Vestibular shield



Patented bilateral blocks



Occlusal plane



U-shaped arch



Key feature

Reduced vestibular shield for optimal comfort with aligners, and a midline upper tab for greater appliance stability. Increased palatal-to-buccal width accommodates the thickness of the aligner.

The lower jaw is positioned edge to edge.

Main actions

- Corrects tongue and lip dysfunctions
- Improves nasal breathing
- Retrains swallowing
- Stabilises occlusal balance

Instructions for use

The bilateral blocks point downwards beneath the tongue. The tongue naturally adopts a raised position.

THE DEVICE MUST BE WORN OVER THE ALIGNER FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

 Prescription and follow-up must be provided by a qualified healthcare professional.



Not suitable for children under 6 years of age or patients with sleep apnoea.

BRACES APPLIANCE

Adjunctive appliance for
fixed brace treatment

Designed to integrate oromyofunctional retraining into orthodontic treatment with fixed braces.



Available in:

M


MEDIUM DENSITY

and

M

MEDIUM DENSITY

WITH BLOCKS



BRACES APPLIANCE

BRACES APPLIANCE

WITH BLOCKS

Also available:

M

MEDIUM DENSITY

UPPER ARCH ONLY

Recommended for tongue retraining when only the upper arch has fixed braces.



Key feature

Upper and lower channels accommodate brackets on the outer surfaces of the teeth.

It supports the action of intermaxillary elastics by positioning the arches edge to edge.

Main actions

- Corrects tongue and lip dysfunctions
- Improves nasal breathing
- Retrains swallowing
- Stabilises occlusal balance
- Relieves temporomandibular joints (TMJs)
- Protects soft tissues

Instructions for use

The bilateral blocks point downwards beneath the tongue. The tongue naturally adopts a raised position.

THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

Prescription and follow-up must be provided by a qualified healthcare professional.




Not suitable for children under 6 years of age, patients with Class III occlusion or prognathic tendencies, or exclusive mouth breathers.

Retention appliances

Maintaining results - Interception

These appliances are used during stabilisation after treatment. They maintain the corrections achieved, prevent relapse and preserve functional and occlusal balance over time. They can also be used during an interceptive phase following an initial stage of treatment.



PROCCLUSEUR

Stabilising results

A retention solution to stabilise the results achieved after interceptive or orthodontic treatment.



How to measure?

Use the dedicated measuring ruler. Measure the upper jaw from the distal surfaces of teeth 22 to 12, at the lower third of the labial surfaces. If the measurement falls between two sizes, choose the larger size.



Learn how to use the measuring ruler and select the most appropriate size for your patient.

Scan the QR code.

Key feature

Indentations from teeth 14 to 24 and 34 to 44. The posterior tracks free molar occlusion to enable physiological function.

Posterior teeth can erupt thanks to the occlusal guidance tracks.

Main actions

- Stabilises results after orthodontic treatment
- Stabilises results after interceptive treatment
- Improves interdigitation of posterior teeth
- Improves nasal breathing
- Helps correct certain deep bites

Instructions for use

Position the Proccluseur in the mouth using the numerical marker on the appliance. The indentations ensure stable, precise and reproducible positioning. Select the size using a dedicated measuring ruler*.

THE DEVICE MUST BE WORN FOR APPROXIMATELY 2 HOURS A DAY AND ALL NIGHT.

 Prescription and follow-up must be provided by a qualified healthcare professional.



Not suitable for children under 36 months or patients with Class III occlusion or sleep apnoea.

* A measuring ruler is available on request when purchasing the Proccluseur.

PROCCLUSEUR

CUSTOMISABLE

Thermoformable appliance
for a personalised fit



How to customise the Proccluseur?



Follow the steps for moulding in the mouth to
achieve a personalised fit and optimal comfort.

Scan the QR code.

Key feature

Creates a customised appliance that accounts for arch width and variations in tooth crown shape and angulation, without laboratory processing. Individualised indentations improve retention in the mouth.

CLINICAL OBSERVATION

At the end of orthodontic treatment, the Proccluseur stabilises occlusion and progressively guides the relationship between the arches. Improved posterior support and better functional integration are observed, reducing the risk of relapse.



DECEMBER 2014



JULY 2015

 Prescription and follow-up must be provided by a qualified healthcare professional.



Not suitable for children under 36 months or patients with Class III occlusion or sleep apnoea.

* A measuring ruler is available on request when purchasing the customisable Proccluseur.

Supporting you beyond the appliances.

Training courses, practical videos, prescribing guidance, clinical support and personalised assistance: Innova Pharm offers a range of resources to simplify your daily practice and support successful treatment outcomes.





Training courses throughout the year.

Who can attend?

Healthcare professionals: orthodontists, general dentists, dental assistants, speech therapists, physiotherapists, osteopaths and allied health professionals.

Our courses are multidisciplinary.



Scan the QR code for
upcoming Innova Pharm
courses, or visit:

www.innovapharm.fr


From
€230 per person*
(lunch and breaks included)

Questions?
We're here
to help.

02 41 58 78 37

mail@innovapharm.fr

www.innovapharm.fr



* For more information about our training courses, visit www.innovapharm.fr.



- ✓ Volume discounts and mixed orders all year (1)
- ✓ 1 FREE appliance for every 10 ordered (2)
- ✓ **Free shipping (3)** on orders of 10+ appliances

For more information,
we're here to help:
mail@innovapharm.fr

 02 41 58 78 37



Innova Pharm
5 bd Pierre Lecoq · BP 30451
49304 CHOLET cedex · France

Innova Pharm 
FUNCTIONAL EXPERTISE SINCE 1999 

Discover all our
appliances online

innovapharm.fr



(1) Mixed orders are available for appliances with the same price. See prices on pages 4-5. (2) Receive one free appliance for every ten ordered. (3) Free shipping within mainland France when ordering 10 or more appliances. Offer applies to appliances only and cannot be combined with other offers.

